

Rb Rabbit mAb


ABZ3417

Product Overview


Name	Rb Rabbit mAb
Accession(Primary)	P06400
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 106 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:50; FC-1:50-1:100; IP-1:20-1:50

Target information (P06400)

Gene Name	RB1
About Gene Symbol	RB1 (RB Transcriptional Corepressor 1) is a protein phosphatase involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Rb in 3T3 lysates using Rb antibody.



Immunocytochemistry analysis of Rb(green) in Hela using Rb antibody, and DAPI(blue)