

TrkB Rabbit mAb


ABZ3418

Product Overview

Name	TrkB Rabbit mAb
Accession(Primary)	P15209
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 92 kDa; Observed MW: 90-140 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:20-1:50; IP-1:20-1:50

Target information (P15209)

Gene Name	Ntrk2
About Gene Symbol	



Western blot analysis of TrkB in mouse brain, rat brain lysates using Trk B antibody.