

SLC31A1 Rabbit mAb

ABZ3419


Product Overview

Name	SLC31A1 Rabbit mAb
Accession(Primary)	O15431
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 21 kDa; Observed MW: 26-34 kDa
Recommended Dilution	WB-1:1000-1:5000


Target information (O15431)

Gene Name SLC31A1

About Gene Symbol



Western blot analysis of SLC31A1/CTR1 in mouse intestines lysates using SLC31A1 antibody.



Western blot analysis of SLC31A1/CTR1 in Hela, CHO-K1, rat Brain, Ramos lysates using SLC31A1/CTR1 antibody.