


Fas (4C3) Mouse mAb

ABZ3420

Product Overview

Name	Fas (4C3) Mouse mAb
Accession(Primary)	P25445
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P25445)

Gene Name FAS**About Gene Symbol**

Western blot analysis of FAS(Cterminus) in Jurkat lysates using FAS(Cterminus) antibody.

