

## 5T4 Rabbit mAb

ABZ3423

### Product Overview

|                             |                    |
|-----------------------------|--------------------|
| <b>Name</b>                 | 5T4 Rabbit mAb     |
| <b>Accession(Primary)</b>   | Q13641             |
| <b>Host Species</b>         | Rabbit             |
| <b>Clonality</b>            | Monoclonal         |
| <b>Isotype</b>              | IgG                |
| <b>Species Reactivity</b>   | Human              |
| <b>Recommended Dilution</b> | IHC-P-1:200-1:1000 |

### Target information (Q13641)

|                          |                                                                                                                                                                                                                                                                                                                                       |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | TPBG                                                                                                                                                                                                                                                                                                                                  |
| <b>About Gene Symbol</b> | TPBG (Trophoblast Glycoprotein) is a protein kinase and signaling regulator involved in Wnt/ $\beta$ -catenin signaling and adaptive immunity. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |