

## ALDH3A1 Rabbit mAb

ABZ3429

### Product Overview

|                             |                    |
|-----------------------------|--------------------|
| <b>Name</b>                 | ALDH3A1 Rabbit mAb |
| <b>Accession(Primary)</b>   | P30838             |
| <b>Host Species</b>         | Rabbit             |
| <b>Clonality</b>            | Monoclonal         |
| <b>Isotype</b>              | IgG                |
| <b>Species Reactivity</b>   | Human              |
| <b>Recommended Dilution</b> | IHC-1:100-1:200    |

### Target information (P30838)

|                          |         |
|--------------------------|---------|
| <b>Gene Name</b>         | ALDH3A1 |
| <b>About Gene Symbol</b> |         |