

AP2 gamma/TFAP2C Rabbit mAb

ABZ3435

Product Overview

Name	AP2 gamma/TFAP2C Rabbit mAb
Accession(Primary)	Q92754
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Recommended Dilution	IHC-P-1:1000-1:2000

Target information (Q92754)

Gene Name	TFAP2C
About Gene Symbol	TFAP2C (Transcription Factor AP-2 Gamma) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.