

ATM Rabbit mAb

ABZ3441

Product Overview

Name	ATM Rabbit mAb
Accession(Primary)	Q13315
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Recommended Dilution	IHC-1:100-1:200

Target information (Q13315)

Gene Name	ATM
About Gene Symbol	ATM (ATM Serine/Threonine Kinase) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.