

Transferrin Receptor Rabbit mAb

ABZ3483

Product Overview

Name	Transferrin Receptor Rabbit mAb
Accession(Primary)	P02786
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Recommended Dilution	IHC-1:100-1:200

Target information (P02786)

Gene Name	TFRC
About Gene Symbol	TFRC (Transferrin Receptor) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.