

## Placental lactogen Rabbit mAb

ABZ3501

### Product Overview

|                             |                               |
|-----------------------------|-------------------------------|
| <b>Name</b>                 | Placental lactogen Rabbit mAb |
| <b>Accession(Primary)</b>   | P0DML3                        |
| <b>Host Species</b>         | Rabbit                        |
| <b>Clonality</b>            | Monoclonal                    |
| <b>Isotype</b>              | IgG                           |
| <b>Species Reactivity</b>   | Human                         |
| <b>Recommended Dilution</b> | IHC-1:100-1:200               |

### Target information (P0DML3)

|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | CSH2 |
| <b>About Gene Symbol</b> |      |