

eIF4E Rabbit mAb

ABZ3515

Product Overview

Name	eIF4E Rabbit mAb
Accession(Primary)	P06730
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Recommended Dilution	IHC-1:100-1:200

Target information (P06730)

Gene Name	EIF4E
About Gene Symbol	EIF4E (Eukaryotic Translation Initiation Factor 4E) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.