

## GDA Rabbit mAb

ABZ3538

### Product Overview

|                             |                 |
|-----------------------------|-----------------|
| <b>Name</b>                 | GDA Rabbit mAb  |
| <b>Accession(Primary)</b>   | Q9Y2T3          |
| <b>Host Species</b>         | Rabbit          |
| <b>Clonality</b>            | Monoclonal      |
| <b>Isotype</b>              | IgG             |
| <b>Species Reactivity</b>   | Human           |
| <b>Recommended Dilution</b> | IHC-1:100-1:200 |

### Target information (Q9Y2T3)

|                          |     |
|--------------------------|-----|
| <b>Gene Name</b>         | GDA |
| <b>About Gene Symbol</b> |     |