

PAX6 Rabbit mAb

ABZ3624

Product Overview

Name	PAX6 Rabbit mAb
Accession(Primary)	P26367
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Recommended Dilution	IHC-1:100-1:200

Target information (P26367)

Gene Name	PAX6
About Gene Symbol	PAX6 (Paired Box 6) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.