

STING Rabbit mAb

ABZ3670

Product Overview

Name	STING Rabbit mAb
Accession(Primary)	Q86WV6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Recommended Dilution	IHC-1:100-1:200

Target information (Q86WV6)

Gene Name	STING1
About Gene Symbol	STING1 (Stimulator Of Interferon Response CGAMP Interactor 1) is an innate immune sensor involved in JAK-STAT/cytokine signaling and innate immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely