

## VIPR1 Rabbit mAb

ABZ3683

### Product Overview

|                             |                  |
|-----------------------------|------------------|
| <b>Name</b>                 | VIPR1 Rabbit mAb |
| <b>Accession(Primary)</b>   | P32241           |
| <b>Host Species</b>         | Rabbit           |
| <b>Clonality</b>            | Monoclonal       |
| <b>Isotype</b>              | IgG              |
| <b>Species Reactivity</b>   | Human            |
| <b>Recommended Dilution</b> | IHC-1:100-1:200  |

### Target information (P32241)

|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | VIPR1 |
| <b>About Gene Symbol</b> |       |