

CD14 Rabbit mAb

ABZ3695

Product Overview

Name	CD14 Rabbit mAb
Accession(Primary)	P08571
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 40 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:500; ICC/IF-1:100-1:200; FC-1:100-1:1000

Target information (P08571)

Gene Name	CD14
About Gene Symbol	