


Caveolin-2 Rabbit mAb

ABZ3698

Product Overview

Name	Caveolin-2 Rabbit mAb
Accession(Primary)	P51636
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:18 kDa;Observed MW: 18 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50

Target information (P51636)

Gene Name CAV2**About Gene Symbol**

Western blot analysis of Caveolin 2 in
HeLa, HUVEC lysates using Caveolin 2 antibody.

