

FAP Rabbit mAb

ABZ3699


Product Overview

Name	FAP Rabbit mAb
Accession(Primary)	Q12884
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:88 kDa;Observed MW: 95 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (Q12884)

Gene Name FAP

About Gene Symbol



Western blot analysis of Fibroblast activation protein, alpha in U-87 MG lysates using Fibroblast activation protein, alpha antibody.

