

CAPON Rabbit mAb


ABZ3700

Product Overview

Name	CAPON Rabbit mAb
Accession(Primary)	O75052
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:56 kDa;Observed MW: 56 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:20-1:50; FC-1:50-1:100

Target information (O75052)

Gene Name	NOS1AP
About Gene Symbol	



Western blot analysis of CAPON in Rat brain
lysates using CAPON antibody.