

HNRNPA0 Rabbit mAb


ABZ3702

Product Overview

Name	HNRNPA0 Rabbit mAb
Accession(Primary)	Q13151
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:31 kDa;Observed MW: 32,34 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:500-1:1000; FC-1:50-1:100

Target information (Q13151)

Gene Name	HNRNPA0
About Gene Symbol	HNRNPA0 (Heterogeneous Nuclear Ribonucleoprotein A0) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of HNRNPA0 in HeLa, A549
lysates using HNRNPA0 antibody.