

MUC1 Rabbit mAb

ABZ3703

Product Overview

Name	MUC1 Rabbit mAb
Accession(Primary)	P15941
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW:122 kDa;Observed MW: 122 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000

Target information (P15941)

Gene Name MUC1**About Gene Symbol**

