

BBS1 Rabbit mAb

ABZ3704


Product Overview

Name	BBS1 Rabbit mAb
Accession(Primary)	Q8NFJ9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:65 kDa;Observed MW: 65 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200

Target information (Q8NFJ9)

Gene Name BBS1

About Gene Symbol



Western blot analysis of BBS1 in SH-SY5Y, 293T lysates using BBS1 antibody.

