

RIP Rabbit mAb


ABZ3705

Product Overview

Name	RIP Rabbit mAb
Accession(Primary)	Q60855
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW:75 kDa;Observed MW: 75 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:100-1:200; FC-1:200-1:500; IP-1:20-1:50

Target information (Q60855)

Gene Name	Ripk1
About Gene Symbol	



Western blot analysis of RIP in RAW
264.7,NIH/3T3 lysates using RIP antibody.