

IL-12A Rabbit mAb


ABZ3706

Product Overview

Name	IL-12A Rabbit mAb
Accession(Primary)	P29459
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW:25 kDa;Observed MW: 35 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; FC-1:100-1:200

Target information (P29459)

Gene Name	IL12A
About Gene Symbol	



Western blot analysis of IL-12A in NIH3T3, Raw 264.7 lysates using IL-12A antibody.