

PRPF8 Rabbit mAb


ABZ3707

Product Overview

Name	PRPF8 Rabbit mAb
Accession(Primary)	Q6P2Q9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:274 kDa;Observed MW: 274 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:2000; ICC/IF-1:50-1:100; FC-1:50-1:100

Target information (Q6P2Q9)

Gene Name	PRPF8
About Gene Symbol	



Western blot analysis of PRPF8 in HeLa, K562
lysates using PRPF8 antibody.