

SHC Rabbit mAb


ABZ3709

Product Overview

Name	SHC Rabbit mAb
Accession(Primary)	P29353
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:63 kDa;Observed MW: 46,52,66 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (P29353)

Gene Name	SHC1
About Gene Symbol	



Western blot analysis of SHC in MCF7, HUVEC
lysates using SHC antibody.