


Angiotensinogen Rabbit mAb

ABZ3710

Product Overview

Name	Angiotensinogen Rabbit mAb
Accession(Primary)	P01019
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:52 kDa;Observed MW: 52 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (P01019)

Gene Name AGT**About Gene Symbol**

Western blot analysis of Angiotensinogen in PC-12 lysates using Angiotensinogen antibody.

