

DTNA Rabbit mAb


ABZ3712

Product Overview

Name	DTNA Rabbit mAb
Accession(Primary)	Q9Y4J8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:84 kDa;Observed MW: 84 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (Q9Y4J8)

Gene Name	DTNA
About Gene Symbol	



Western blot analysis of DTN-A in Human fetal
lysates using DTN-A antibody.