

PPAR gamma Rabbit mAb


ABZ3713

Product Overview

Name	PPAR gamma Rabbit mAb
Accession(Primary)	P37231
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW:58 kDa;Observed MW: 58 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100-1:200; FC-1:200-1:1000

Target information (P37231)

Gene Name	PPARG
About Gene Symbol	



Western blot analysis of PPAR gamma in Human fetal lysates using PPAR gamma antibody.