

AP2S1 Rabbit mAb


ABZ3714

Product Overview

Name	AP2S1 Rabbit mAb
Accession(Primary)	P53680
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:17 kDa;Observed MW: 17 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; IP-1:10-1:100

Target information (P53680)

Gene Name	AP2S1
About Gene Symbol	AP2S1 (Adaptor Related Protein Complex 2 Subunit Sigma 1) is a signaling adaptor protein involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate



Western blot analysis of AP2S1 in SH-SY5Y ,293T
lysates using AP2S1 antibody.