

## C11B2/CYP11B2 Rabbit mAb

ABZ3717

### Product Overview

|                             |                                          |
|-----------------------------|------------------------------------------|
| <b>Name</b>                 | C11B2/CYP11B2 Rabbit mAb                 |
| <b>Accession(Primary)</b>   | P19099                                   |
| <b>Host Species</b>         | Rabbit                                   |
| <b>Clonality</b>            | Monoclonal                               |
| <b>Isotype</b>              | IgG                                      |
| <b>Species Reactivity</b>   | Human                                    |
| <b>Calculated MW</b>        | Calculated MW:58 kDa;Observed MW: 45 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IP-1:10-1:100          |

### Target information (P19099)

**Gene Name** CYP11B2**About Gene Symbol**

