

EHMT2 Rabbit mAb


ABZ3719

Product Overview

Name	EHMT2 Rabbit mAb
Accession(Primary)	Q96KQ7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:132 kDa;Observed MW: 160-180 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:500-1:1000; FC-1:50-1:100; IP-1:20-1:50

Target information (Q96KQ7)

Gene Name	EHMT2
About Gene Symbol	EHMT2 (Euchromatic Histone Lysine Methyltransferase 2) is a histone methyltransferase involved in chromatin remodeling and epigenetics. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to



Western blot analysis of EHMT2/G9A in HEK-293, Jurkat lysates using EHMT2/G9A antibody.