

Nucleolin Rabbit mAb


ABZ3720

Product Overview

Name	Nucleolin Rabbit mAb
Accession(Primary)	P19338
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:77 kDa;Observed MW: 100 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200; FC-1:200-1:2000

Target information (P19338)

Gene Name	NCL
About Gene Symbol	



Western blot analysis of Nucleolin in Human placenta lysates using Nucleolin antibody.