

## ADAR1 Rabbit mAb


ABZ3721

### Product Overview

|                             |                                                      |
|-----------------------------|------------------------------------------------------|
| <b>Name</b>                 | ADAR1 Rabbit mAb                                     |
| <b>Accession(Primary)</b>   | P55265                                               |
| <b>Host Species</b>         | Rabbit                                               |
| <b>Clonality</b>            | Monoclonal                                           |
| <b>Isotype</b>              | IgG                                                  |
| <b>Species Reactivity</b>   | Human                                                |
| <b>Calculated MW</b>        | Calculated MW:136 kDa;Observed MW: 150 kDa           |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:50-1:100;<br>FC-1:10-1:100 |

### Target information (P55265)

|                          |                                                                                                                                                                                                                                                                                                                                      |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | ADAR                                                                                                                                                                                                                                                                                                                                 |
| <b>About Gene Symbol</b> | ADAR (Adenosine Deaminase RNA Specific) is a cytokine and immune signaling protein involved in JAK-STAT/cytokine signaling and innate immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of ADAR1 in HeLa lysates  
using ADAR1 antibody.