

MKK7 Rabbit mAb

ABZ3722

Product Overview

Name	MKK7 Rabbit mAb
Accession(Primary)	O14733
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:48 kDa;Observed MW: 48 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; FC-1:50-1:100

Target information (O14733)

Gene Name	MAP2K7
About Gene Symbol	

