


PGC1 beta Rabbit mAb

ABZ3723

Product Overview

Name	PGC1 beta Rabbit mAb
Accession(Primary)	Q86YN6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:113 kDa;Observed MW: 113 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (Q86YN6)

Gene Name PPARGC1B**About Gene Symbol**

Western blot analysis of PGC1 beta in HeLa, MCF-7 lysates using PGC1 beta antibody.

