

GAA Rabbit mAb


ABZ3724

Product Overview

Name	GAA Rabbit mAb
Accession(Primary)	P10253
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:105 kDa;Observed MW: 76 ,95,110kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P10253)

Gene Name	GAA
About Gene Symbol	



Western blot analysis of GAA in 293T, Human placenta lysates using GAA antibody.