

STAT3 Rabbit mAb


ABZ3725

Product Overview

Name	STAT3 Rabbit mAb
Accession(Primary)	P40763
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:88 kDa;Observed MW: 88 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100; FC-1:200-1:500

Target information (P40763)

Gene Name	STAT3
About Gene Symbol	STAT3 (Signal Transducer And Activator Of Transcription 3) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of STAT3 in HeLa, SH-SY5Y
lysates using STAT3 antibody.