


PPIC Rabbit mAb

ABZ3726

Product Overview

Name	PPIC Rabbit mAb
Accession(Primary)	P45877
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW:23 kDa;Observed MW: 23 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:100-1:200

Target information (P45877)

Gene Name PPIC**About Gene Symbol**

Western blot analysis of PPIC in HCT116 lysates using PPIC antibody.

