


## ADRA1B Rabbit mAb

ABZ3727

### Product Overview

|                             |                                          |
|-----------------------------|------------------------------------------|
| <b>Name</b>                 | ADRA1B Rabbit mAb                        |
| <b>Accession(Primary)</b>   | P35368                                   |
| <b>Host Species</b>         | Rabbit                                   |
| <b>Clonality</b>            | Monoclonal                               |
| <b>Isotype</b>              | IgG                                      |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                        |
| <b>Calculated MW</b>        | Calculated MW:57 kDa;Observed MW: 57 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:500      |

### Target information (P35368)

**Gene Name** ADRA1B**About Gene Symbol**

Western blot analysis of ADRA1B in HeLa, PC-3 lysates using ADRA1B antibody.

