


CAD Rabbit mAb

ABZ3728

Product Overview

Name	CAD Rabbit mAb
Accession(Primary)	P27708
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW:243 kDa;Observed MW: 243 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:200-1:500

Target information (P27708)

Gene Name CAD**About Gene Symbol**

Western blot analysis of CAD in 293 lysates using CAD antibody.

