

Gephyrin Rabbit mAb


ABZ3730

Product Overview

Name	Gephyrin Rabbit mAb
Accession(Primary)	Q9NQX3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:80 kDa;Observed MW: 93 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (Q9NQX3)

Gene Name	GPHN
About Gene Symbol	GPHN (Gephyrin) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Gephyrin in 293T lysates using Gephyrin antibody.