

SOX1 Rabbit mAb

ABZ3734

Product Overview

Name	SOX1 Rabbit mAb
Accession(Primary)	O00570
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:39 kDa;Observed MW: 39 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100; FC-1:100

Target information (O00570)

Gene Name SOX1**About Gene Symbol**

