


KLF10 Rabbit mAb

ABZ3735

Product Overview

Name	KLF10 Rabbit mAb
Accession(Primary)	Q13118
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:53 kDa;Observed MW: 53 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q13118)

Gene Name KLF10**About Gene Symbol**

Western blot analysis of KLF10 in Jurkat, K652 lysates using KLF10 antibody.

