


DARS Rabbit mAb

ABZ3736

Product Overview

Name	DARS Rabbit mAb
Accession(Primary)	P14868
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:57 kDa;Observed MW: 57 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:20-1:50

Target information (P14868)

Gene Name DARS**About Gene Symbol**

Western blot analysis of DARS in Jurkat, K562
lysates using DARS antibody.

