


C1orf93 Rabbit mAb

ABZ3737

Product Overview

Name	C1orf93 Rabbit mAb
Accession(Primary)	Q8TBF2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:21 kDa;Observed MW: 21 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (Q8TBF2)

Gene Name PRXL2B**About Gene Symbol**

Western blot analysis of C1orf93 in MCF7 lysates using C1orf93 antibody.

