


Syntaxin 3 Rabbit mAb

ABZ3738

Product Overview

Name	Syntaxin 3 Rabbit mAb
Accession(Primary)	Q13277
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:33 kDa;Observed MW: 33 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (Q13277)

Gene Name STX3**About Gene Symbol**

Western blot analysis of Syntaxin 3 in HepG2, Raw 264.7 lysates using Syntaxin 3 antibody.

