


INTS10 Rabbit mAb

ABZ3739

Product Overview

Name	INTS10 Rabbit mAb
Accession(Primary)	Q9NVR2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:82 kDa;Observed MW: 82 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q9NVR2)

Gene Name INTS10**About Gene Symbol**

Western blot analysis of INTS10 in HeLa, Jurkat lysates using INTS10 antibody.

