


NIPP1 Rabbit mAb

ABZ3740

Product Overview

Name	NIPP1 Rabbit mAb
Accession(Primary)	Q12972
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:39 kDa;Observed MW: 39 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q12972)

Gene Name PPP1R8**About Gene Symbol**

Western blot analysis of NIPP1 in HeLa, Jurkat lysates using NIPP1 antibody.

