


MRPS35 Rabbit mAb

ABZ3742

Product Overview

Name	MRPS35 Rabbit mAb
Accession(Primary)	P82673
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW:37 kDa;Observed MW: 37 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P82673)

Gene Name MRPS35**About Gene Symbol**

Western blot analysis of MRPS35 in HeLa ,Jurkat
lysates using MRPS35 antibody.

